

# KL-602 Syringe Pump

## Service Manual



Beijing KellyMed Co., Ltd.

## Importance

This manual provides information needed to service the Syringe Pump. This manual is intended for use by certified biomedical engineering technicians or personnel with equivalent experience and training in servicing this type of equipment. It is recommended that the user complete the training class geared specifically to the Syringe Pump.

## **Revision History**

The following revision history table summarizes changes in this document.

Revision No.	<b>Revision Date</b>	Description of Changes	
1.0	21/05/2010/	Initial version.	
1.1	16/04/2015	Update pictures of front cover, front view and rear view.	
1.2	21/12/2021	English grammar correction.	
1.2	31/12/2021	Update EC representative of back cover.	

## Contents

1. Description of pump	1
1.1 Specifications	
1.2 Front view	
1.3 Rear view	
1.4 Syringe view	
2. System block diagram	
3. Tools and equipments (for	service and test)5
4. Troubleshooting	
5. Check power supply	
6. Replace battery	
7. Replace FPC cable betwee	n motor board & inter-connection board15
8. Adjust the syringe detector	board
9. Check software version	
10. Syringe calibration	
11. Occlusion calibration	
12. Software upgrade	
12.1 Tools and program.	
12.2 Connect PC and put	mp with RS23226
12.3 Software upgrade	

## 1. Description of pump

#### 1.1 Specifications

POWER REQUIREMENTS	230V AC, 50Hz 8 * AA NIMH battery 9.6VDC		
BATTERY LIFE	1.5~2 years		
POWER CONSUMPTION	$\leq 15 \text{VA}$		
FLOW RATE	0.1 ~ 99.9mL/h       (in 0.1mL/h increments)         100 ~ 1300mL/h (in 1mL/h increments)         0.1~400mL/h for 10mL syringes.         0.1~600mL/h for 20mL syringes.         0.1~800mL/h for 30mL syringes.         0.1~1300mL/h for 50mL syringes.		
ACCURACY	±2%		
VOLUME	0.001 ~ 99.999ml (in 0.001ml increments) 100.00 ~ 999.99ml (in 0.01ml increments) 1000.0 ~ 9999.9ml (in 0.1ml increments)		
DELIVERY LIMIT	0.1 ~ 99.9ml (in 0.1ml increments) 100 ~ 9999ml (in 1ml increments)		
KVO	0.1mL/h~1.0ml/h (user selected)		
SYRINGE SPECIFICATION	10ml, 20ml, 30ml, 50ml		
ANTI-BOLUS	Release pressure when pump get occluded		
BOLUS	1 ml ~ 5ml		
BOLUS RATE	Maximum Rate of each syringe type		
OPERATION CONDITIONS	5~40°C, 20~90% RH (no condensation)		
STORAGE CONDITIONS	-20~55°C, ≤93% RH (no condensation)		
WARRANTY	1 years		

Table 1.1 Specifications

#### 1.2 Front view



Figure 1.1 Front view

1	LCD DISPLAY	2	CHARGE LED	3	AC LED
4	ALARM LED		WORKING LED	6	[ON/OFF]KEY
7	[START/STOP] KEY		[SILENCE]KEY	9	[MENU]KEY
10	[MODE]KEY	11	[PURGE] KEY	12	[UP]KEY
13	[DOWN]KEY	14	[ENTER]KEY	15	[SHUTTLE]KEY
16	SYRINGE INSTALL SLOT	17	SYRINGE PLATE	18	WATERPROOF CLOTH
19	PRESSURE SENSOR	20	OPEN NUT	21	PLUNGER HOOK
22	PUSHER	23	V SHAPE SLOT		

#### 1.3 Rear view



Figure 1.2 Rear view

1	CASCADE SCREWS	2	CASCADE BAR	3	PUMP FEET
4	BATTERY COVER	5	MODEL TAG	6	PUMP HANDLE
7	COMMUNICATE PORT	8	AC POWER INLET	9	CLAMP SCREW
10	CLAMP FOLDING SWITCH	11	POLE CLAMP	12	HAND WHEEL

1.4 Syringe view



Figure 1.3 Syringe view



## 2. System block diagram

Figure 2.1 System block diagram

1	Power	Power cord
2	Replace Battery	New battery (8 * AA NIMH battery 9.6VDC) Cross-screwdriver
3	Replace fuses	Cross-screwdriver An electric soldering iron One fuse (2A 250V)
4	Disassembly	Cross-screwdriver

## 3. Tools and equipments (for service and test)

## 4. Troubleshooting

No.	The symptom of failure	Cause	Solution	
		Alarm when the pump working	<ol> <li>Calibrate pressure sensor (See Chapter 11)</li> <li>Check and see if there is stuff in the pipe(remove)</li> </ol>	
1	Occlusion alarm		3、Check the plane of PRESSURE SENSOR(Front view 1.1 (19)) Whether above shell(if not, replace Pressure Sensor)	
			4. Check with the value of Pressure sensor change when press it with finger(if not, replace Pressure Sensor)	
	The set value is too sensitive		<ol> <li>Adjust the sensitive of occlusion alarm</li> <li>Calibrate pressure sensor(See Chapter 11)</li> </ol>	
2 Occlus	Occlusion but not alarm	Bad connection in the circuit board	Check the connecting point and reinstall in the original status(see Chapter 7)	
		Bad Pressure sensor	Replace a new one	
		Failure of the system	Replace a new main board	
3 Abnormal LCD & L display	Abnormal LCD & LED display	Bad LCD	Replace a new one	
		The indicator fails to work	Replace a new lamp or send back the circuit board and repair	
4	Key fails to work	Bad key	Replace a new key	

_		The battery out	Connect to the AC and charge the battery
5	Low battery alarm	No AC display or low battery alarm	Check the power supply circuit(see Chapter 5)
		Motor is weak or fail to work	Replace a new motor
6	The pump fail to work	Failure in the Driver IC or circuit board	Check the circuit board and replace a new driver IC
7	Injection accuracy big error	Injection accuracy above 2%	Repeat defining a custom syringe.(see Chapter 10)
8		Install syringe error	repeat install syringe Press the [PURGE] key twice, and the second keep pressing is about two seconds
0		Bad Pressure sensor	Replace a new one
0	Nut amon alarma	NUT no installation position	Press the [PURGE] key twice, and the second keep pressing is about two seconds until the "Nut Err" disappeared on LCD
9	Nut error alarm	Bad connection in the circuit board	Check the connecting point and reinstall in the original status(see Chapter 7)
10	0Syringe PlateSyringe detect board positionDetection errorerror		Adjust the syringe detect board position(see Chapter 8)
		AC power cord is not inserted properly	Check the AC power cord connection.
11	Pump cannot be turned on.	Internal battery is out	Stop the operation of the pump and replace with a new battery through your local authorized dealer.
		The voltage of the internal	Recharge the battery fully for more than 8 hours by connecting the pump to an AC power
		battery is low.	outlet and turn it on.

## 5. Check power supply

Connect AC plot

- Step 1 If 230VAC absent, Please check the power supply socket connection to the board.
- Step 2 If not 230V, please replace fuse.
- Step 3 If not 14.75V-15V, please replace the power supply board.
- Step 4 If not above 9.0V, please charging for the internal battery.
- Step 5 If not about 5.0V, please replace the main board.





#### Other boards in pump



Jog& shuttle board



Motor driver board



Syringe detect board

## 6. Replace battery

Step 1 New battery (Internal battery must be supplied by the manufacturer)



Step 2 Open the cover of battery



Step 3 Pull off the conjunction, take the failed battery out



Step 4 Replace a new battery in the battery position



Step 5 Put the conjunction into the slot



Step 6 Close the cover and fix it.

## 7. Replace FPC cable between motor board & inter-connection board

Step 1 Prepare a new FPC cable (Must be supplied by manufacturer)Step 2 Open the shell



Step 3 Take off the defective FPC cable



Step 4 Replace a new FPC cable





Step 5 Close the shell and fix it

## 8. Adjust the syringe detector board

Step 1 Open the shell

Step 2 Scroll the screws off with a cross screwdriver



Step 3 Adjust the position of syringe detector board in or out

Step 4 Test the detector function status

Turn on the Syringe Pump then press [START/STOP]KEY into the bellowed interface



In this interface loaded different syringe size (10ml 20ml 30ml 50ml), can display syringe size correctly, if not display correctly, return STEP3, Adjust the position of syringe detector board in or

out slightly Until finding the correct position.

Step 5 After the detector function adjustment, fix the screws on.

Step 6 Close the shell and fix it

#### NOTE:

If adjust the position of syringe detector board in or out failed, please enlarge the holes in which the screw in by filing then copy the above steps.

## 9. Check software version

Through the following steps to check software version

- 1) Turn on the pump
- 2) Press [MENU] KEY to enter the main menu

Rate Mode(ml/h)	
VOL-Time Mode	
Weight Mode	▼
MENU	BAT.

3) Select "Special Setting" from main menu and press "enter" key to enter "Special Setting" menu

Default Setting		
Log		
Product Info.	$\checkmark$	
Special	BAT.	

4) Select "Product Info" from "Special Setting" menu and press "enter" key to enter Product Information Display



## 10. Syringe calibration

The Syringe Pump support custom two brand : custom 1, custom 2, Each of four types: syringe including 10ml, 20ml 30ml 50ml

Tools:

- 1. Vernier caliper, or ruler
- 2. syringe

Ready:

1) Measure A (barrel length) and B (plunger length) and record.



- 2) Install syringe (Example 10 ml)in pump (According to the instructions of the installation method of installation).
- 3) Into custom syringe interface.
- 4) When the brand for black when we show, Press [ENTER]KEY. Let custom 1 or 2 against the white custom, the diagram below:

Brand:	CUSTOM 1		
A 76.0	0mm	В	28.00mm
Confirm			Cancel
Custom			BAT.

5) Choose the name of self-define syringe(self-define 1 and self-define 2) by turning shuttle key. This selection custom 2, the diagram below:

Brand	: CUSTO	DM 2	
A 76	.00mm	В	28.00mm
Confirm			Cancel
Custom			BAT.

6) Press [ENTER]KEY. After confirm, turning shuttle key or use [UP]KEY [DOWN]KEY choose A . show as the diagram below:

Bra	nd: CUST	OM 1	
A	76.00mm	В	28.00mm
Cor	nfirm		Cancel
Custom			BAT.

7) Input the value of A and Press [ENTER]KEY confirm.

Such as measuring the length of a syringe effective scale for 75.8 mm, syringe push handle length b for 25.8 mm A study is adjusted for 75.8, adjust the numerical values of item B for 25.8

- 8) Adopt the same method, set the data item B.
- 9) Confirm input. Syringe custom finish.
- 10) Into syringes options, check the definition of syringe parameters are now stored. As shown in figure

Brand:	CUSTOM2		
10ml	75.80	25.80	
20ml	00.00	00.00	
30ml	00.00	00.00	
50ml	00.00	00.00	
syringe		BAT.	

If you want to continue to define different types of syringe repeat the step 1) - step 8).

### 11. Occlusion calibration

Through the following steps to calibrate the pressure sensor:

- 1) Turn on the pump.
- 2) Press [MENU] KEY to enter the main menu.



3) Select "Special Setting" from main menu and press "enter" key to enter Special Setting menu.



4) Select "Maintenance" from Special Setting menu and press "enter" key to enter system Maintenance Interface.(For example calibrate pressure sensor for 10ml).



5) Install syringe then connect syringe and manometer



6) Press [START/STOP] KEY of the pump when it displays as above picture shows. then motor of the pump starts running. meanwhile "PUMP-1" changes to "PUMP-2".

AC: 367	BAT: 3	75
P:342	PUMP-2	S10
K: 20	MARK1	SAVE
Maintain.		BAT.

7) Press [START/STOP] KEY of the pump When manometer displays 0.02MPa ,then press [DOWN] KEY Select "MARK1"as shown below.

AC: 367 BAT: 375				
P:342	PUMP-1	S10		
K: 20	MARK1 SAVE			
Maintair	۱.	BAT.		

8) Press [START/STOP] KEY of the pump when display above. then motor of the pump starts running. While "MARK1" into "MARK2".

K: 20	MARK2	SAVE	
AC: 36	7 BAT: 37	′5	
P:342	PUMP-1	S10	

 Press [START/STOP] KEY of the pump When manometer displays 0.10MPa ,then press [DOWN] KEY Select "SAVE" as shown below.

Maintain.	BAT.
K: 20 MARK1	SAVE
P:342 PUMP-1	S10
AC: 367 BAT: 375	5

10) Press [ENTER] KEY of the pump save the record. Calibrate pressure sensor for 10ml syringe finished.

Repeat the above step4 to step10 to calibration pressure sensor for 20ml, 30ml and 50ml syringe

## 12. Software upgrade

- 12.1 Tools and program
- 1) RS232 Serial Port Connector Line.



2) Problem for Software upgrade



3) Software (new version)



- 12.2 Connect PC and pump with RS232
- 2.1) Operating system required
- Win 9X/Me /2000/XP/ (Windows XP professional edition is recommended)
- 2.2) Connect PC and the Syringe Pump with RS232.



## WARNING

- A Disconnect the AC power cord.
- B Turn off the Syringe Pump before connection

#### NOTE:

A If there is not RS232 communication port, replace it by a "USB/RS232 line". Connection diagram as below.



B USB to serial line.



12.3 Software upgrade

Step 2

Step 1 Connect PC and Syringe Pump with RS232.

M	BootLoader Control. exe
F L	[TUDU: \File desc.

Double-click this icon

, Open the software upgrade

tool--"Boot Loader Control.exe", The tool's interface is as below.

AN2295 BootLoader fo:	The file path of softwar	e new version
	1	V Open S19 file
Communication port(COM1)		Rescan
Communication Port	COM Speed: 9600	~
semiconductor	BaudRate	otload Exit

Step 3 Select COM Speed

Select the "COM Speed" value "115200".as below.

🏯 AN2295 BootLoader fo	HC08/HCS08	Control - F	···· 🔳 🗖 🔀
			V Open S19 file
Communication port(COM1)			Rescan
	COM Speed:	Bootload	Exit

Step 4 Select the Software program file (New version).

Click the "Open S19 file" button , Select the file name "L602.01.05.abs.s19".

🛃 AN2295 Bo	otLoader for	r HCO8/HCSO8 Co	ontrol - Fi		
打开				V Open 519 file	?
查找范围(I):	🗁 Bootloader	_GUI	<b>v</b> (	3 🜶 🖻 🛄 -	
1000 我最近的文档	EL602.01.05.	abs. s19			
<b>[]</b> 桌面					
》 我的文档					
<b>夏</b> 夏 我的电脑					
國金融 國金融 國金融 國金融 國金融 國金融 國金融 國金融 國金融 國金融 國金融 國金融 國金融 					
	文件名 (J): 文件类型 (J):	L602.01.05.abs.s1 S19 files (*.s19;	9 *. sx;*. srec)	*	打开 (1) 取消

Step 5 The interface of selecting the software program file is as below.



Step 6 Click the button "Bootload", Enter downloading interface.



Step 7 Press the [ON/OFF] KEY of the pump (Keep holding, don't release) until waiting process is finished. Software upgrade is process as below.



Please release the [ON/OFF] KEY of the pump when the interface closed automatically, software upgrading is finished.

NOTE:

2)

1) If your computer has a firewall, please remove the restrictions on the upgrade tool



when you use "USB-RS232 line" to

replace RS232, Open the software upgrade tool--"Boot Loader Control.exe", The tool's interface is shown as below.

🛔 AN2295 BootLoader for HC08/HC	508 Control - Fr 🔳 🗖 🔀
E:\KL-602\Documents\Bootloader_GUI\L602.01.05	.abs.s19 🗸 Open S 19 file
USB-SERIAL (COMn)	Rescan
COM Spee	d: 115200 💌
"n "may ha 1 to 2	Bootload Exit

3) If failed to upgrade procedures, ( Such as: Release the [ON/OFF]KEY of the pump in the

process of upgrading or RS232 etc. ), Please repeat the above step 1 to step 7 to upgrade again.

4) If upgrade procedures is failed the pump will not startup.

#### NOTE:

Occlusion alarm occura frequently after the software upgrade, please follow the method to calibrate pressure sensor.



Beijing KellyMed Co., Ltd. Room 115, 4 # Building, No. 2 Tiyuan West Road, Haidian District, Beijing, 100084, China Tel: +86-10-82490385 Fax: +86-10-65587908

**Umedwings Netherlands B.V.** 



Treubstraat 1, 2288EG, Rijswijk, The Netherlands Tel: +31(0)642758955 E-mail: ar@umedwings.eu