



# SK-ET75077A

Emergency Trolley



 **Medspares**  
Pacific  
EQUIPMENT  
SPARES  
ACCESSORIES

Medspares Pacific Ltd  
PO Box 2070 Raumati Beach 5032  
Wellington - New Zealand  
P: +64 4 299 7610 - F: +64 4 299 2784  
e: sales@medspares.co.nz - www.medspares.co.nz

**SAIKANG** / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE  
Stock Code: 870098

## QUALITY CONTROL

### MATERIAL

Steel and aluminum structure inside, ABS and antibacterial materials outside, easy to clean.



### DRAWER

Drawer insize big & deeper!  
Inner with partitions can be well organized

### EXTENDABLE WORKTOP

Laterally extendable worktop



### TABLE-BOARD

With handles on both sides, surface with transparent soft glass

### IV POLE

Movable IV pole can exchange position with defibrillator shelf



### DEFIBRILLATOR SHELF

The plate surface is made of virgin ABS material. The position of the IV pole and the defibrillator shelf can be exchanged for easy use.

### POWER OUTLET

The power socket is movable and detachable, Easy to use, Located above the CPR board to convenient the use of oxygen cylinders



### OXYGEN BOTTLE HOLDER

Inserted oxygen bottle holder, belts fixed on trolley main body

### TRASH CAN

Two ABS trash cans with lid



### WHEELS

TPR tire no worn out after running 30KM, save anti-winding hard shell, United forming without bolts  
Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500times

### DEBRIS BASKET

With two debris baskets on both sides. The big one is for two sharp containers and the small one is for sanitizer

