

M3A Vital Signs Monitor



Portable Design

Infrared Ear Thermometer



Medspares Pacific
EQUIPMENT
SPARES
ACCESSORIES

PO Box 2070 Raumati Beach 5032
Wellington - New Zealand

P: +64 4 299 7610 - F: +64 4 299 2784

e: sales@medspares.co.nz - www.medspares.co.nz

EDANUSA

www.edanusa.com

M3A Vital Signs Monitor

EDAN M3A Vital Signs Monitor provides essential measurement of patient's SpO₂, NIBP and Temperature. It offers an effective solution and brings high performance vital signs monitoring at your finger tips.

30,000 NIBP Measurement Review	72 h Trend Review	800 Alarm Review
--	-----------------------------	----------------------------

- 3.5" high resolution color TFT-LCD screen & LED display
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 17 hours of continuous work
- Compatible with Patient Care Viewer PC data management software
- Nurse call
- Suitable for adult, pediatric and neonatal patients
- Barcode scanner support (optional)
- Ethernet connection for bi-directional communication with MFM-CMS central monitoring system (optional)
- Wi-Fi support with optional AP



Mounting Solutions



Wall Mount

Rolling Stand

Bed Rail Clamp

Specification

Configuration	SpO ₂	NIBP	Quick Temp	Infrared Ear Thermometer
M3A(SpO ₂ only)	✓	—	—	—
M3A(NIBP only)	—	✓	—	—
M3A(SpO ₂ +NIBP)	✓	✓	—	—
M3A(NIBP+Quick Temp)	—	✓	✓	—
M3A(NIBP+Infrared Ear Thermometer)	—	✓	—	✓
M3A(SpO ₂ +NIBP+Quick Temp)	✓	✓	✓	—
M3A(SpO ₂ +NIBP+Infrared Ear Thermometer)	✓	✓	—	✓



4901 Morena Blvd., Suite 505 | San Diego, CA 92117 | Tel: 888.850.4597