## NEW VITALOGRAPH micro

The Ultimate Low Cost, High Performance Spirometer





Data you can rely on. People you can trust.


## NEW VITALOGRAPH micro


The Ultimate Low Cost, High Performance Spirometer

The Vitalograph<sup>®</sup> micro<sup>™</sup> is a high quality hand-held spirometer at a great value price from the leading provider of spirometry devices to general healthcare and occupational health markets. This reliable, easy-to-use, entry-level model offers all required functionalities in a full feature spirometer. The micro is unrivalled in its price range and is the ideal choice for fast, accurate testing in environments where mobility is important.

The micro features the same high performance measuring technology used across the Vitalograph range and requires no user calibration. In contrast to turbine spirometers the micro does not feature any sensitive moving parts, guaranteeing robustness and the highest possible measuring accuracy - even at very low or high flow rates.

- Low cost, entry-level model from the leader in general healthcare and occupational health spirometry test devices
- Durable measurement technology
- No need for expensive disposable sensors, turbines, spirettes or flow tubes
- No user calibration required
- Optimal hygiene via Vitalograph Bacterial Viral Filters (BVF™)
- Touchscreen for simple and intuitive operation
- Includes Vitalograph Reports Software for creating editable PDF reports (for printing, filing and uploading to electronic medical record (EMR) systems)
- Features latest GLI spirometry predicted equations and Z-scores





✓Low Running Costs 
✓Optimal Hygiene

Ordering Information

Authorised Distributor

## Accessories

63310 Vitalograph micro

28350 Vitalograph BVF Bacterial Viral Filters (50)20303 Noseclips (200)

Vitalograph Ltd, UK. Maids Moreton, Buckingham, MK18 1SW, England. T: (01280) 827110 F: (01280) 823302 E: sales@vitalograph.co.uk www.vitalograph.co.uk Vitalograph Ltd, Intl. Maids Moreton, Buckingham, MK18 1SW, England. T: +44 1280 827120 F: +44 1280 823302 E: sales@vitalograph.co.uk www.vitalograph.eu Vitalograph Inc, 13310 West 99th Street Lenexa, Kansas, 66215, USA. Toll Free: 800 255 6626 T: (913) 730 3200 F: (913) 730 3232 E: vitcs@vitalograph.com www.vitalograph.com

Vitalograph®, micro™, and BVF™ are trademarks or registered trademarks of Vitalograph Ltd.