



iM3
Vital Signs Monitor

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,
Kengzi Subdistrict, Pingshan District, Shenzhen | 518122, P.R. China
+86.755.26898326 | www.edan.com.cn | info@edan.com.cn

© 2017 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-iM3-V1.0-20170928

A Slim Monitoring Experience

iM3 Vital Signs Monitor



Smarter, Larger, yet Slimmer

Combining a stunning ultra-slim design, a large capacitive touch screen, together with a rich set of leading functions, iM3 takes the lead and presents itself as the beginning of a new generation of vital signs monitors. Bringing in full options of SpO₂, NIBP and Quick Temp/Infrared Ear Temp and differentiated work modes, iM3 serves for both long time monitoring and spot check applications.



Touch Screen



- **Standard Parameters**
EDAN SpO₂, EDAN NIBP, PR
- **Optional Parameters**
Nellcor™ SpO₂, SunTech® NIBP, EDAN Quick TEMP, Infrared Ear TEMP

Redefinition of Adaptive Work Modes

- **Round Mode**
Special work mode designed for daily ward rounding application. Worklist could be easily input through USB disk or downloaded from EMR system. Moreover, the color-coded patient list comes with square windows interface design, you can tell patients' status and information in a quick glance.



Patient management



Patient status with color code

- **Monitoring Mode**
For continuous patient monitoring with real-time data, real-time alarms, and history trend.
- **Spot Mode**
For spot-check applications, such as out-patient, emergency, etc. Real-time vitals are acquired and displayed, no alarm will be seen or heard.

TEMP Module Options



Infrared Ear TEMP



EDAN Quick TEMP



Rolling Stand

Non-volatile built-in Memory

Spot & Round Mode

16,000,000
Measurements

Monitoring Mode (for each patient file)

200	1200	120 h
Alarm	NIBP Measurements	Trend