


Your Vision Is Our Goal

# SonoEye Vet



## SonoEye Vet Ultrasound Device

---

Fast, flexible and easy to assist  
diagnosis at the point of care.

# SonoEye Vet handheld ultrasound

☑ Lightweight | ☑ Waterproof | ☑ Thumb-play | ☑ Portable


**This cost effective solution, with no additional costs, gives you access to advanced ultrasound techniques at the point of care.**

The adaptable, portable, design works on any device; whether a tablet or mobile phone, making it possible for veterinarians to support animal patient to patient measurements any time, anywhere. SonoEye provides excellent scanning experience for animals of different species.

- **Half the weight of a modern smartphone, SonoEye features a compact, highly integrated design that is easy to hold.**
- **SonoEye is waterproof and can be safely cleaned and sterilised to assist infection control.**
- **Cost effective solution with no additional costs.**
- **With simplified workflows SonoEye is easy to use on your smartphone or tablet; either IOS or Android. One touch can start and finish an ultrasound exam.**
- **SonoEye, integrated with the latest technology, can provide multiple solutions for the scanning of different species:**
  - **Feline**
  - **Canine**
  - **Ovine**
  - **Equine**
  - **Other Species**



# QUALIFIED CLINICAL PERFORMANCE



**Canine Cardiac, C Mode**



**Canine Bladder, B Mode**



**Canine Aorta, C Mode**



**Canine Kidney, B Mode**



**Canine Liver B Mode**



**Canine Aorta, B Mode**

# SonoEye offers a variety of veterinarian solutions with its specialist scanners, helping to meet your every need.



## SonoEye V2

7.5Mhz, 40mm, Linear  
Abdomen, Vascular,  
Intestine, Bladder,  
Adrenal Gland, etc.



## SonoEye V3

3.0Mhz, Phased Array  
Cardiac, Abdomen,  
Reproduction, etc.



## SonoEye V5

3.5Mhz, Convex  
Abdomen,  
Reproduction, etc.



## SonoEye V6

6.0Mhz, Micro Convex  
Abdomen, Vascular,  
Bladder, Intestine,  
Adrenal Gland, Cardiac  
EM, etc.



SonoEye Vet can help to provide high quality user experience in any location.



Gain



Depth



TGC



## Multiple accessories available



Arm bag



Belt bag



Messenger bag



Carry Case



Stand



Sterile probe cover and gel



**CHISON**  
Value Beyond Imaging

CHISON Medical Technologies Co. Ltd  
No.9, Xinhuihuan Road, Xinwu District, Wuxi,  
Jiangsu, China 214028.



**Medspares  
Pacific**  
EQUIPMENT  
SPARES  
ACCESSORIES

**Medspares Pacific Ltd**  
PO Box 2070 Raumati Beach 5032  
zWellington - New Zealand  
P: +64 4 299 7610 - F: +64 4 299 2784  
e: sales@medspares.co.nz - www.medspares.co.nz