



# Hugemed VLM(VLRM/VLDM)

Video laryngoscope



# HugeMed VLRM/VLDM

When time is critical, make your go-to device the HugeMed VLRM/VLDM.

Ideal for demanding OR, ED, and ICU settings, the VLRM/VLDM provides a consistently clear view of the airway, and enhanced maneuverability enabling fast intubations. The VLRM/VLDM system gives you a choice of 13 reusable blades and 5 disposable blades to cover a wide range of patients—from infant to morbidly obese.

Add real-time recording and an onboard tutorial, and you'll see why HugeMed VLRM/VLDM is the one you go to for advanced video laryngoscopy.

## VLRM



## Reusable Blades

### VLRM/VLDM=ADVANCED AIRWAY VIEWS

- Digital color monitor with LCD-quality resolution
- Reusable video laryngoscope with digital-resolution camera
- Innovative blade sizing, angulation, and camera positioning
- Anti-fog feature, with a rapid heating profile to resist lens clouding and secretions

## VLDM



## Disposable Blades

### VLRM/VLDM=ADVANCED USER FEATURES

- Integrated, real-time recording helps confirm tube placement and facilitate teaching
- Ergonomic design and rugged construction for multiple clinical applications
- Reusable blades in 13 sizes:  
MILLER 00, 0, 1,  
MAC 1,2,3,4,5  
MACD 1,2,3,4,5
- Also available in 5 single-use sizes

## Advanced Video Laryngoscopy:

Real-time recording. LCD clarity.



### Applications

- First-use intubations, replacing direct laryngoscopy
- Normal or restricted oropharyngeal views
- Cormack–Lehane Grades I–IV laryngeal views
- Trauma airways—excellent when dealing with blood and secretions in the airway
- Video-guided tube exchange in the ICU

### Specifications

Type	VLRM	VLDM
Display	8 Inch LCD monitor HD Display	
Resolution	1024*768RGB	
Data audiovisual output	Support	
Camera Resolution	2M Pixel	
View Angle	66°	
Illumination	≥800LUX	
Battery Life Time	120–150min	
Battery Capacity	10000mAh	
Voltage	3.7V	
Blade type	Reusable Blades	Disposable Blades
Material of Blade	SUS316L	PC
Shot and Record Function	Support	
Anti-fog	Anti-fog when it's on	
Blade Sizing	Miller 00/0/1; Mac 1/2/3/4/5; MacD 1/2/3/4/5	Mac 1/2/3/4/5
HDMI&USB Port	Include	
Connector for Screen and Handle	Include	
Bracket (with tray)	Include	
Attachments	Adapter, Power cord, Data connection line	
Blades Disinfection	Refer to the reusable laryngoscope instructions	Refer to the Disposable laryngoscope instructions



### Silicone Cover Blades

Mac 1/2/3/4/5

Option: blades with channel (Mac 3/ 4)



### Reusable Blades

Model No.		Size(L*W*H)		Applicable patients
Miller 00		72X14.4X29.5		premature
Miller 0		84X15X29.5		newborn
Miller 1		110X15X29.5		infant
Mac 1	MacD 1	101X18X34	96X18X34	infant
Mac 2	MacD 2	115X21X36.5	115X22X50	children
Mac 3	MacD 3	130X24X40	127X22X34	adult
Mac 4	MacD 4	136X23X43	136X22X57	difficulty (obese)
Mac 5	MacD 5	150X23X51	147X22X59	difficulty



### Disposable Blades

Model No.	Size(L*W*H)	Applicable patients
Mac 1	95X28X12.5	infant
Mac 2	110X29X13	children
Mac 3	120X27X13	adult
Mac 4	135X22X14	difficulty (obese)
Mac 5	120X29X16	difficulty

# ACCESSORIES



## MOBILE STAND

A stand-alone configuration for use in a variety of hospital settings



## RIGID STYLET

For use with endotracheal tubes 6,0 mm and larger.

Recommended to facilitate successful intubations.



## EMS BAG

Color: BLUE

Size: 8inch



Shenzhen HugeMed Medical Technical Development Co.,Ltd

Add: 8F,Block A, Tempus Building,1st Qingshuihe Road, Luohu District

Shenzhen City, Guangdong, China



**Medspares Pacific Ltd**

PO Box 2070 Raumati Beach 5032

zWellington - New Zealand

**P:** +64 4 299 7610 - **F:** +64 4 299 2784

**e:** [sales@medspares.co.nz](mailto:sales@medspares.co.nz) - [www.medspares.co.nz](http://www.medspares.co.nz)