

EndoView BR Series Single use Flexible Bronchoscope



Technical Specifications

| Model | | BR-M22 | BR-M32 | BR-M40 | BR-M52 | BR-M50 | BR-M58 |
|-----------------|---------------------------|----------|--------|--------|--------|--------|--------|
| Optical system | Field view | 120° | | | | | |
| | Depth of field | 3-50mm | | | | | |
| Insertion tube | Distal end outer diameter | 2.2mm | 3.2mm | 4.0mm | 5.2mm | 5.0mm | 5.8mm |
| | Working length | 600mm | | | | | |
| Working Channel | | / | 1.2mm | 2.0mm | 2.6mm | 2.8mm | 3.4mm |
| Angulation | Bending | 210° U/D | | | | | |
| Total Length | | 855mm | | | | | |

| Display | |
|------------|-----------------|
| Screen | LCD |
| | Size: 8 inch |
| Function | Photo |
| | Video |
| | White Balance |
| | Freeze |
| Battery | Type: Lithium |
| | Maximum 4 hours |
| Connection | HDMI |

| Single use Flexible Bronchoscope | |
|----------------------------------|---------|
| Order Information | Product |
| 2016B-PA005 | BR-M22 |
| 2016B-PA004 | BR-M32 |
| 2016B-PA003 | BR-M40 |
| 2016B-PA001 | BR-M50 |
| 2016B-PA002 | BR-M52 |
| 2016B-PA006 | BR-M58 |

Shenzhen HugeMed Medical Technical Development Co.,Ltd

Add: 416-1,516-1, Building 2, No. 1, Mawu Road, Baoan Community, Yuanshan Street, Longgang District, Shenzhen 518115, China
 Tel: +86-755-2227 5833
 Fax: +86-755-2227 5833
 Email: sale@hugemed.net

www.hugemed.net

AN ENDOSCOPE SOLUTION COMPANY


EndoView BR Series
Single use Flexible Bronchoscope

- Always Available
- Always Brand New
- No Costly Repair
- No Contamination



EndoView BR Series Single use Flexible Bronchoscope

- Plug in and ready to go
- Multiple size to satisfy various clinical settings
- Navigate with ease




High Quality Image



210° / 210° high deflection,
Ensure full navigation

- Slim Model: OD 2.2/ ID 0 mm
OD 3.2/ ID 1.2 mm
- Routine Model: OD 4.0/ ID 2.0 mm
OD 5.2/ ID 2.6 mm
- Therapy Model: OD 5.0/ ID 2.8 mm
OD 5.8/ ID 3.4 mm



Customized Two buttons:
WB / Photo / Video / Freeze

USB and HDMI connection,
easy data transfer



Compatible with conventional
accessories