

# Automatic Sealing Machine SEAL1000

- 7 Inches Big Touch Screen User Interface
- Rapid Temperature Heat Up & Cool Down
- Easy Operation & Automatic Continuous Sealing
- Meet The Demand Of Big Clinic
- Customized Setup



- Sealing speed: 10± 0.5m /min
- Sealing margin: 0 ~ 35 mm adjustable
- Sealing line width: > 12 mm
- Operating temperature: 60 ~ 220°C adjustable
- Temperature control accuracy: ≤ 1%
- Printing method: needle printing
- Working environment: 10 ~ 40°C
- Ac power supply: 220V±10% 50Hz
- Power: 500 WMaximum
- current: 3.2a
- Fuse : 5A×2
- Overall dimensions: 560×260×220mm

