

A98-4 ELECTRIC OBSTETRIC BED



SAIKANG/Medical

Healthier your life

A98-4

FEATURES AND HIGHLIGHTS

Height adjustment







Technical parameters

 External Size (LxW) 	1980x960mm(±30mm)
 Height 	650-940mm(±20mm)
■ Back-rest Adjustment	0-60°
 Reverse Trendelenburg 	0-7°
 Leg-rest Adjustment 	-9°-35°, leg section is removable
Weight capacity	200KG



Technical configuration

✓ 5" Double side castors, central locking✓ Central locking pedal2pcs
☑ Central locking pedal 2pcs
☑ ABS head board with safe lock 1pc
✓ ABS side rail with angle indicator 1set
✓ Hand grip 1pair
☑ Leg holder 1pair
☑ Basin 1pc
☑ Hand controller 1pc
☑ Power line 1pc
✓ Leg strip 1 pair
☑ IV pole prevision 2pcs







Material

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean.



Bed surface

Bed surface by plastic steel, made by virgin material. The cover is detachable, and easy to clean.



Head Board and Guardrail

The Head board and side rail made by virgin ABS material. The head board is movable.



Hand controller

Equipped with one hand controller.



Handle

The hand bar can be dismounted and movable.



Lithotomy rod

The leg support can be dismounted and movable.



Caster

5 inch central controlled caster for moving obstetric table.
TPR tire, quiet and wear-resistant,

TPR tire, quiet and wear-resistant, can be washed and sterilized by high temperature, save anti-winding hard shell, united forming without bolts.



Accessory

Equipped with stainless steel sewage basin and other various accessories.



Mattress

The mattress adopts high density memory sponge and is once formed. Mattress is movable.



Foot section

Foot section can be retracted, convenient for gynecological and delivery surgical.



















Design

Manufacturing

Maintenance



South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina TEL: +55 (11) 2758-0438 Rua Universal, Hollywood, São Bernardo do Campo, Brazil

TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower, upper hill, nairobi, Kenya TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd

Sala Thammasop Subdistrict Thawi Watthana District Bangkok Thailand TEL: +66 642741046

UK Branch

Uk Healthward International Limited

DEPT 111, 196 High Road Wood Green, London, N22 8HH United Kingdom TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG' 454004, Russian Federation, Chelyabinsk,

Universitetskaya Naberezhnaya, Russia TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC

Musaffah P.O Box 53200, Abu Dhabi, UAE Tel: +971 585710652



PO Box 2070 Raumati Beach 5032 zWellington - New Zealand P: +64 4 299 7610 - F: +64 4 299 2784 e: sales@medspares.co.nz - www.medspares.co.nz