

SK-AT75077A2

Anesthesia Trolley



LEADING BRAND OF CHINA MEDICAL FURNITURE Stock Code: 870098

QUALITY CONTROL



DRAWER

Drawer insize big & deeper! Centralized lock, Inner with partitions can be well organized

EXTENDABLE WORKTOP

Laterally extendable worktop





With a sharp container at one side

SHARP CONTAINER

TRASH CAN Two different color ABS trash cans with one lid, the lid controlled by foot pedal with damping function









ANESTHESIA BOX With five pcs Anesthesia box each line, In total two lines



WHEELS

TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts Pass dynamic test: Bearing 120kg run 30KM and pass obstacles 500times



