



SK-CT75077C2

Treatment Trolley



Medspares Pacific Ltd
PO Box 2070 Raumati Beach 5032
Wellington - New Zealand
P: +64 4 299 7610 - F: +64 4 299 2784
e: sales@medspares.co.nz - www.medspares.co.nz

SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

QUALITY CONTROL

TABLE-BOARD

ABS upper top with extended up thin rims on three sides, Surface with transparent soft glass



MATERIAL

Steel and aluminum structure inside, ABS and antibacterial materials outside, easy to clean.



DRAWER

With ABS inserts inside, which can be separated freely as required.



DEBRIS BASKET

With three debris baskets on both sides. One side is for two sharp containers and another side is for sanitizer



WHEELS

TPR tire no worn out after running 30KM, Save anti-winding hard shell, United forming without bolts
Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500time



TRASH CAN

Two ABS trash cans with lid

