

# X36 EXAMINATION TABLE

Registration code: SKD-A



SAIKANG/Medical

Healthier your life

# X36

# FEATURES AND HIGHLIGHTS

# Technical parameters

External Size (LxWxH)	1860x690x460-910 mm
Back Plate Upward Fold	≥80°
Back Plate Downward Fold	≥50°
<ul> <li>Leg Plate Upward Fold</li> </ul>	≥65°
<ul> <li>Leg Plate Downward Fold</li> </ul>	≥45°
Bearing Capacity	≥150kg
<ul> <li>Headrest Hole</li> </ul>	150 mm
<ul> <li>Input Voltage</li> </ul>	110-230v
Output Power	24W

# Technical configuration

✓ Motor	1pc
<b>☑</b> Pillow	1pc
☑ Gas Spring	2pcs
☑ Brake Pedal	4pcs
☑ 3" Double Wheel	4pcs
✓ Manual Operator	1pc

# Optional

☐ Bed Surface Foot Lift	1pc
☐ Roll Up Tissue Holder	1pc















#### **Bed Frame**

The bed frame surface is painted with 1.5mm powder coated steel, and welded by Panasonic Robotic seamlessly. Strong corrosion resistance, environmentally-friendly, bright surface, maintain an overall beautiful appearance. With this structure, the whole bed maintains more stable and makes patients feel more safe.



#### Mattress

The bed platform is integrally made of the high-end composite leather and hand-made cotton mattress, of which the curve design is ergonomic, elastic and soft, breathable, environmentally-friendly, easy to be cleaned, anti-bacterial, impactresistant



#### **Pneumatic Pole Support**

Two air springs are equipped on the back and legs to provide users with a quiet and safe environment during diagnosis and avoid interference.



#### **Control System**

CE certified linear actuator motor, one-buttom free lift, low noise, and long service life. Manual (foot control) electric overall lifting, manual back-rest lifting, leg-rest lifting adjustment, height adjustment.



### **Hand Controller**

The bed is equipped with Hand-held controller, so that the optimal height is easy to be adjusted by pressing the single button gently. The hook on the back is convenient to take down for patients and nurses.



## Bed foot equipment

Equipped with 4pcs 3-inch double-sided wheels, 2 fuctions side brake with 4 pedals, which are convenient to move and fix. The design of 360° rotating pulley made of polyurethane PU is strong, durable, nonsip, strong in pressure resistance, and good in stability.



### **Assorted Pillow**

The neoteric design of the pillow is narrow on the top and wide on the bottom. Its well-shaped and elegant design matches with modern aesthetic ideas and can blend into the environment harmoniously. Pillow can be disassembled according to demand.



#### **Round Covered Edge**

The bed platform was covered in 360°. The soft-touch of surface would be safer for users.



## Design

Bed surface two fold three ergonomic design, adjust the height according to the patient's height, personalized design help diagnosis. Dual use of bed and chair, high bearing lasting stability, no fear of heavy pressure impact, to ensure safe use.

















Design

Research

Manufacturing

Maintenance

Service



# South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina TEL: +55 (11) 2758-0438 Rua Universal, Hollywood, São Bernardo do Campo, Brazil

TEL: +55 (11) 2758-0438

# East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower, upper hill, nairobi, Kenya TEL: +254714084806

#### South East Asia Branch

Saikang Medical Technology Ltd

Sala Thammasop Subdistrict Thawi Watthana District Bangkok Thailand TEL: +66 642741046

## UK Branch

Uk Healthward International Limited

DEPT 111, 196 High Road Wood Green, London, N22 8HH United Kingdom TEL: +31 (0) 202111106

# Russia Branch

LLC 'SAIKANG'

454004, Russian Federation, Chelyabinsk, Universitetskaya Naberezhnaya, Russia TEL: +7 (919)300-00-06

# Mid East Branch

Saikang Medical Supplies (Mid East) LLC

Musaffah P.O Box 53200, Abu Dhabi, UAE Tel: +971 585710652





Medspares Pacific Ltd
PO Box 2070 Raumati Beach 5032
Wellington - New Zealand
P: +64 4 299 7610 - F: +64 4 299 2784
e: sales@medspares.co.nz - www.medspares.co.nz